



## **COLLECTION FLAWS**

### Fatal Flaws

- Specimen ID numbers do not match on CCF and bottle.
- Both collector signature and printed name are missing.
- Bottle seal is broken, absent, or shows evidence of tampering.
- Neither primary, nor split specimen, have at least 30 mL of urine.
- CCF and specimen did not arrive in the same package.

### Correctable Flaws

- Just collector's signature is omitted.
- No donor signature and no 'Remarks'.
- Donor's printed name is missing.
- Just collector's printed name is missing.
- Certifying scientist's signature is missing on a non-negative.
- Non-federal form used in error for DOT test.
- Expired version of CCF used.
- Drug tests to be performed not checked.
- Date of collection missing.
- Specimen temperature was not checked, and there are no 'Remarks'.

### Minor Flaws

- Address or phone number errors.
- Reason to test is not checked.
- Donor ID information incorrect or omitted.
- Testing authority box not checked.
- Split specimen was not checked.